## Table of Contents

<b>Thermosaver Hypothermic Rescue Bag Processing Instruction -ENG-.....</b>	<b><a href="#">2</a></b>
Read before use / Important information.....	2
Operational lifespan of Thermosaver Hypothermic Rescue Bag.....	2
Processing between uses.....	3
Identification of consumable components.....	3
Processing instruction .....	4
Warnings and limitations.....	10
<b>Thermosaver Hypothermic Rescue Bag - Instruksjoner For Vask -NOR- .....</b>	<b><a href="#">13</a></b>
<b>Thermosaver Hypothermic Rescue Bag Instruktioner för tvätt -SWE- .....</b>	<b><a href="#">24</a></b>
<b>Thermosaver Hypothermic Rescue Bag -Brugsanvisning -DAN- .....</b>	<b><a href="#">35</a></b>

# Thermosaver Hypothermic Rescue Bag Processing Instruction - ENG -

Manufacturer: Thermosaver AS

Device(s): TSR-5431R, TSR-5431MR, TSR-5431MJKR,  
TSR-5531, TSR-5531M, TSR-5531MJK

## Read before use / Important information

**Personnel responsible for the processing of the Thermosaver Hypothermic Rescue Bag, as assigned by the organization, shall be made aware of the content of this instruction, and should be followed based on the organizations internal processes and procedures, prior to first use and first processing.**

The Instructions for Use of the device, provided in physical copy, shall always be followed and read before use, in conjunction with this instruction before processing, available in downloadable PDF format at the manufacturer's website. Internet connection and a modern web browser is required capable of opening PDF files.

Thermosaver Hypothermic Rescue bag, (ref. Device(s) model numbers) is a **multiple patient multiple use** medical device intended by the manufacturer to be reused on multiple patients for multiple uses. Based on the Spaulding classification it is considered a non-critical reusable medical device, based on the established low degree of risk for infection involved in relation to the intended use of the device(s).

The device(s) is not placed on the market under sterile packaging conditions and is not intended for sterile use, neither is intended to be sterilized before use.

The device(s) requires processing between patients and uses following the recommendations within this instruction. The maximum number of allowable uses for the Hypothermic Rescue Bag is **9 processing cycles**.

## Operational lifespan of Thermosaver Hypothermic Rescue Bag

The Thermosaver has a maximum lifetime of 5 years or 10 times use (first use + 9 reuses), whichever comes first. An extension of the lifespan after solely storage or minimal usage cycles is the responsibility of the manufacturer.

Damaged, torn, degraded device(s) resulting from unforeseen circumstances shall never be used and represents a risk for patients. If you are uncertain about the condition of your device, do not use it before consulting the manufacturer for confirmation on its safe use.

## Processing between uses


Before the second use of the Thermosaver Hypothermic Rescue bag the steps detailed within this processing instructions must be followed.

This instruction is intended to be followed by healthcare professionals.

## Identification of consumable components

The “Vapour barrier and absorbent” is single use only and **shall be discarded after each use**, following environmentally sustainable practices.

Compatible consumables listed within this section are available for purchase from the distributor.

### Plastic Bag:

### Compression Bag (recommended to be used for storage after processing):

\*Colour may vary based on model

### Hood:

\*Colour may vary based on model

## Processing instruction

The steps covered within this section represent a complete processing cycle and covers:

- Preparation before processing
- Cleaning, disinfection and drying
- Inspection and maintenance
- Preparation for use and how to replace consumable components
- Packaging and storage

### Preparation before processing

Follow these steps to prepare the Hypothermic Rescue Bag for safe handling and/or further processing:

- 1) A used Thermosaver Hypothermic Rescue Bag should be kept in the plastic bag provided in the pocket of the bag until the start of the cleaning procedure. The “**Vapour barrier and absorbent**” should be disposed of as hospital waste directly after use. **If the bag is wet, the cleaning procedure needs to be started within 24 hrs of use.**
- 2) Before cleaning, disinfection and drying, the Thermosaver Hypothermic Rescue Bag should be examined:
  - For damages: A damaged bag should be disposed of.
  - For soiling: Remove general soiling such as dirt and debris from the bag.
  - For severe soiling: Severe soiling would normally be limited to the hood. A severely soiled hood should be disposed of and replaced.
- 3) Verify that the number processing cycles already performed has not exceeded the maximum number of 9 allowed cycles. If the Thermosaver exceeds maximum processes, it should be returned to the manufacturer for disposal and recycling.
- 4) Check pockets and remove any objects from the pockets, including the IFU before storing the bag for cleaning, disinfection and drying.

### Cleaning, disinfection and drying

Cleaning and disinfection of the Thermosaver must follow the guidelines applicable to hospital fabrics washable at 60 °C including disinfection according to national guidelines at an approved laundry facility.

The cleaning, disinfection and drying process must be conducted by a laundry that holds the necessary approvals for cleaning of hospital linen including, but not limited to, certification of the disinfection process. The laundry shall be free to follow their internal established and validated processes in compliance with national regulations, for defining the concentration of disinfectant and detergent, the mix ratio and temperatures, duration, exposure time as part of the washing cycle, with considerations to choose a method suitable for hospital fabrics washable at 60 °C.

- The Thermosaver bag can be tumble dried at 70 °C to 80 °C.
- **DO NOT IRON.**


- The use of an automated cleaning, disinfection and drying process is recommended which can be achieved by the use of an industrial washing machine (subject to validation requirements by the laundry), suitable for this purpose.

## Inspection and Maintenance

### Inspection

The Thermosaver Hypothermic Rescue bag must be inspected after cleaning, disinfection and drying and prepared for use. If any damage to the bag is observed, or any of the below inspection points does not pass, the bag should be returned to the manufacturer.

To identify parts referenced within the inspection check points refer to the below diagrams:

Inspection check points	Suitable for use	Not suitable for use
<p>Check that all the Velcro's (red and black) are clean of debris.</p>	 <p>No debris, foreign objects are observed on the Velcro's surface</p>	 <p>The surface of the Velcro's is contaminated with debris.</p>

Inspection check points	Suitable for use	Not suitable for use
<p>Check the Inner skin, and Outer skin of the fabric material to be free of ruptures</p>	 <p>No visible tearing, ruptures are noticed.</p>	 <p>Tearing or ruptures are observed on the fabric.</p>
<p>Check all the carrying handles seems and the handle's material to be free of tears or ruptures</p>	 <p>No visible tearing or ruptures were noticed, neither on the handles material or on the handles seems.</p>	 <p>One or more handles or seems present signs of tearing or rupture.</p>
<p>Check the content of the label on the device to be complete and legible with special attention for the serial number.</p>	 <p>All information on the label is present, legible, and readable, including the serial number.</p>	 <p>Information on the label is washed out, parts of it is not legible or readable anymore.</p>
<p>Serial label remains attached to the bag</p>	 <p>The label is attached to the bag. Minor degradation of attachment on the corners is still acceptable.</p>	 <p>Complete detachment of the label</p>

## Maintenance

After successful inspection after the processing a record must be kept, to track the number of processing cycles performed on the serial number of the Thermosaver Hypothermic Rescue Bag. The organization owning the bag is responsible for maintaining a track record by serial number. It is required that the record is stored for safekeeping.

You can locate the serial number (SN) on the label, according to the picture.

The serial number for traceability is in the following format depending on the model:

TSR-5431/XXXXX

TSR5431MJK/XXXXX

where XXXXX is a 5-digit number representing your bags unique number.

Thermosaver provides as part of this instruction a template that shall be filled in with the below details and shall be traceable to the serial number of the processed Hypothermic Rescue Bag:

- What is the serial number of the processed Hypothermic Rescue Bag?
- When the Thermosaver bag was cleaned, disinfected and dried?
- Which approved laundry performed the cleaning, disinfection and drying?
- Which vapour barrier was fitted?
- Which packaging method was used (vacuum packaging vs compression bag)?

## Preparation for use and how to replace consumable component

Prior to packaging, the Thermosaver Hypothermic Rescue Bag must be equipped with a new Vapour Barrier and a hood if it has been disposed of as part of the Preparation before processing.

### 1) Fitting of the Vapour Barrier:

- Lay down the Hypothermic Rescue Bag, fully opened on a clean, flat surface
- Position the vapour barrier on the inside area covering entirely the shape of the unit with the number decals facing upwards.

Additional resource in the form of an instructional video (first 10 seconds) showcasing the placement of the vapour barrier is available online <https://thermosaver.no/> under Instructional Videos section. Requires internet connection and modern web browser.

### 2) Fitting of the Hood

- On the top (wider section) of the bag attach the new hood to the Velcro.

## Packaging and Storage

### Packaging using the compression bag

- Ensure that the vapour barrier is fitted before packaging.
- Log roll the bag from the sides. Have the bag laid down on a clean, flat surface for this operation.

- Start folding it from the leg section, 3-4 times, until a compact, compressed size is achieved.

- Place the folded bag within the compression bag.

Additional resource in the form of an instructional video showcasing the packaging within the compression bag is available online <https://thermosaver.no/> under Instructional Videos section. Requires internet connection and modern web browser.

### Storage

Store the processed Thermosaver Hypothermic Rescue bags in dry, clean, and dust-free environment.

If the bag cannot be stored under dry conditions, it is recommended that the processed bag is vacuum packaged for longer term storage.

**Warning: Storage of the processed bag in a humid environment, may affect the water absorption capabilities of the vapour barrier!**

## Warnings and limitations

The below information represents the process details, which has been tested by Thermosaver to determine the maximum number of allowable reuses. This process is provided as recommendation to be followed by operators engaging in the processing of Thermosaver Hypothermic Rescue bag.

All warning, limitations, precautions listed within these instructions shall be carefully read and followed.

<b>Limitations on processing</b>	It is the obligation of the user to ensure that the processing procedure, including resources, materials and personnel, is suitable to achieve the required results. State-of-the-art technology and national laws require compliance with validated processes.
<b>Information on validation of processing</b>	Thermosaver relied upon on validated processes provided by the laundry Sentralvaskeriet Østlandet AS where the washing for cleaning, disinfection tests has been performed with the following process and equipment details:  Cleaning agents: Cool Chemistry from Christeys Dosage: as per laundry's internal certification; Washing Equipment: Girbau HS4012 washer-extractor Wash Program: with pre-wash, main-wash and rinse cycle; Water Temperature: 60 °C Wash cycle duration: 60 minutes Drying method and temperature: Tumble Drying at 70 °C to 80 °C max. Drying Duration: 30 minutes. G Factor: 135 G
<b>Warnings</b>	If the above-mentioned cleaning agents and machines are not available, or the cleaning, disinfection and drying cycle is different, the user is obliged to validate his procedure accordingly.

## Disposal

Thermosaver can be contacted for the correct disposal of the bag at the end of its shelf life or after the expiration of its in-use life. The Thermosaver Hypothermic Rescue Bag is to be completely recycled by the manufacturer.

## **Manufacturer contact information**

Thermosaver AS  
Rosenholmveien 25,  
1414 Trollåsen, Norway

Telephone: +47 928 92 186

Mail: [post@thermosaver.no](mailto:post@thermosaver.no)

Web: [www.thermosaver.no](http://www.thermosaver.no)

Monday – Friday      8 – 16

Saturday – Sunday    Closed

## Thermosaver Hypothermic Rescue Bag Processing Tracker

**Serial number:**

<b>Processing Count</b>	<b>Date</b>	<b>Approved Laundry Facility performing cleaning, disinfection and drying</b>	<b>Person name performing Inspection</b>	<b>Packaging method chosen – vacuum vs. compression</b>	<b>Which Vapour barrier was fitted</b>
1					
2					
3					
4					
5					
6					
7					
8					
9					

# Thermosaver Hypothermic Rescue Bag - Instruksjoner For Vask -NOR-

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## Drifts levetid for Thermosaver Hypothermic Rescue Bag

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